 Sontec Instruments®

SPINE LINE

**Black McCullough Retractor Set**


**Black McCullough Retractor Set**


- SI-8001 set
- 416-150 retractor w/extended rack 8"
- SI-575-800 m-black case
- SI-2.04 narrow blades 4cm (pair)
- SI-2.05 narrow blades 5cm (pair)
- SI-2.06 narrow blades 6cm (pair)
- SI-2.07 narrow blades 7cm (pair)
- SI-2.74 wide blades 4cm (pair)
- SI-2.75 wide blades 5cm (pair)
- SI-2.76 wide blades 6cm (pair)
- SI-2.77 wide blades 7cm (pair)
- SI-H3 hook 3cm
- SI-H4 hook 4cm
- SI-H5 hook 5cm
- SI-H6 hook 6cm
- SI-MH3 multi-tooth hook 3cm
- SI-MH4 multi-tooth hook 4cm
- SI-MH5 multi-tooth hook 5cm
- SI-MH6 multi-tooth hook 6cm



416-150







Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS

Sontec Instruments, Inc.  
7248 South Tucson Way  
Centennial, Colorado 80112 USA  
Toll-free / 800.821.7496  
Local / 303.790.9411  
Fax / 303.792.2606

[SontecInstruments.com](http://SontecInstruments.com)