



 Sontec Instruments®

KERRISON RONGEUR

Universal Klick Kerrison Rongeur Set with Tray

404-608-SET Set

Classic Standard Set includes:

404-6081	1 mm	40° 8" shaft
404-6082	2 mm	40° 8" shaft
404-6083	3 mm	40° 8" shaft
404-6084	4 mm	40° 8" shaft
404-6085	5 mm	40° 8" shaft
404-6086	2 universal handles	

18-3564 tray only

To configure you own personalized set, refer to "Kerrison Section" > "Kerrison Universal Rotating 360° Rongeur" for all possible options

*different size shafts available up request



Special

Kerrison Universal Rotating 360° Rongeur

404-606B handle *special order

Bayonet shaft / black ceramic / 12 1/2" (31.8 cm)

404-6061B 1 mm bite 40°
 404-6062B 2 mm bite 40°
 404-6063B 3 mm bite 40°

Endoscopic / black ceramic / 16" (40.6 cm)

404-6063EB 3 mm bite 40°

Quick trigger release / shaft locks in 8 positions through a 360° rotation / 8" (20cm) shaft

404-6082B 2 mm bite 40° up curved with tissue ejector
 404-6082C 2 mm bite 40° down curved with tissue ejector
 404-60825E 2.5 mm bite 40°
 404-6083B 3 mm bite 40° up curved with tissue ejector
 404-6083E 3 mm bite 40° endoscopic



404-6083B



404-6083E



404-60825E



404-6082B



404-6082C




404-6062B



404-6063B



404-6063EB



Handle sold
separately
from shaft

360 degree rotation

404-606B

Kerrison Universal Rotating 360° Rongeur

Quick trigger release / with tissue ejector / shaft locks in 8 positions through a 360° rotation / 8" (20cm) shaft / 40°

404-606	handle
404-605	full set includes: handle all 40° sizes 8", all 90° sizes 8"
404-6081	1 mm bite
404-6082	2 mm bite
404-6083	3 mm bite
404-6084	4 mm bite
404-6085	5 mm bite

With tissue ejector / 10" (25.5cm) shaft

404-606	handle
404-607	full set- includes: handle all 40° sizes 10", all 90° sizes 10"
404-6086	1 mm bite
404-6087	2 mm bite
404-6088	3 mm bite
404-6089	4 mm bite
404-6090	5 mm bite

7" (17.5cm) shaft

404-6091	0.8 mm bite
404-6093	1.5 mm bite



404-6081



404-6082



404-6083



404-6084



404-6085



404-6091




404-6093



with tissue ejector

404-6082



404-6082

360 degree rotation

404-606

Handle sold
separately
from shaft

Kerrison Universal Rotating 360° Rongeur

Quick trigger release / with tissue ejector / shaft locks in 8 positions through a 360° rotation / 8" (20cm) shaft / 90°

404-606	handle
404-605	full set includes: handle all 40° sizes 8", all 90° sizes 8"
404-6071	1 mm bite 90°
404-6072	2 mm bite 90°
404-6073	3 mm bite 90°
404-6074	4 mm bite 90°
404-6075	5 mm bite 90°

With tissue ejector / 10" (25.5cm) shaft

404-606	handle
404-607	full set- includes: handle all 40° sizes 10", all 90° sizes 10"
404-6076	1 mm bite 90°
404-6077	2 mm bite 90°
404-6078	3 mm bite 90°
404-6079	4 mm bite 90°
404-6080	5 mm bite 90°



404-6071



404-6072



404-6073



404-6074




404-6075



with tissue ejector

404-6073



404-6073

360 degree rotation

Handle sold
separately
from shaft

404-606

Kerrison Universal Rotating 360° Rongeur Handle

404-606 stainless steel
404-606B black ceramic
404-606SB black ceramic silicone handle



Sontec Beatty Foraminotomy Rongeur

7 3/4" (19.5 cm) shaft / curved up

404-461 2 mm bite

404-462 3 mm bite

404-463 4 mm bite

7 3/4" (19.5 cm) shaft / curved down

404-450 2 mm bite


404-451 3 mm bite

404-452 4 mm bite


* Also available with rounded footplate



slightly curved up



curved up



curved down



Beatty Foraminotomy Rongeur

**Black ceramic coating / 8" (20.3 cm) shaft /
slightly curved up**

404-461B 2 mm bite
404-462B 3 mm bite
404-463B 4 mm bite

* Also available with rounded footplate



Kerrison 40° Rongeur


Black ceramic coated / thin
footplate punch / 10 mm opening

		Bite Width	Length of Shaft		
			6" / 15 cm	7" / 18 cm	8" / 21 cm
40° Up-Angle		1 mm	404-490B	404-539B	404-400B
		1.5 mm		NEW SIZE 404-5391	
		2 mm	404-491B	404-540B	404-402B
		2.5 mm		NEW SIZE 404-5401	
		3 mm	404-492B	404-541B	404-404B
		4 mm	404-493B	404-542B	404-406B
		5 mm	404-494B	404-543B	404-408B
		6 mm	404-495B	404-544B	404-410B



Kerrison 90° Rongeur

Black ceramic coated / 10 mm opening

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
 Up-Bite	1 mm	404-470B	404-519B	404-412B
	2 mm	404-471B	404-520B	404-414B
	3 mm	404-472B	404-521B	404-416B
	4 mm	404-473B	404-522B	404-418B
	5 mm	404-474B	404-523B	404-420B
	6 mm	404-475B	404-524B	404-422B



404-520B

Kerrison 90° Rongeur

Thin footplate-punch / ceramic coating


	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
90° Down-Bite 	1 mm	404-480B	404-530B	404-424B
	2 mm	404-481B	404-525B	404-426B
	3 mm	404-482B	404-526B	404-428B
	4 mm	404-483B	404-527B	404-430B
	5 mm	404-484B	404-528B	404-432B
	6 mm	404-485B	404-529B	404-434B



Kerrison 40° Laminectomy/Lumbar Rongeur

Black ceramic coating / 8" (20 cm) shaft / wide mouth 15 mm jaw opening / with & without ejection pin

	Bite Width	w/Ejector Pin	wo/Ejector Pin	wo/Ejector Pin Stainless Steel
40° Up-Angle 	1 mm		404-337B	404-332
	2 mm	404-332B	404-338B	404-333
	3 mm	404-333B	404-339B	404-334
	4 mm	404-334B	404-340B	404-335
	5 mm	404-335B	404-341B	



Kerrison 40° Rongeur

Thin footplate-punch / 10 mm opening

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
40° Up-Angle 	1 mm	404-490	404-539	404-400
	2 mm	404-491	404-540	404-402
	3 mm	404-492	404-541	404-404
	4 mm	404-493	404-542	404-406
	5 mm	404-494	404-543	404-408
	6 mm	404-495	404-544	404-410



404-402

Kerrison 90° Rongeur

Thin footplate punch / 10 mm opening / up

		Bite Width	Length of Shaft		
			6" / 15 cm	7" / 18 cm	8" / 21 cm
Up-Bite		1 mm	404-470	404-519	404-412
		2 mm	404-471	404-520	404-414
		3 mm	404-472	404-521	404-416
		4 mm	404-473	404-522	404-418
		5 mm	404-474	404-523	404-420
		6 mm	404-475	404-524	404-422



404-520

Sontec Kerrison Rongeur 90° Down Bite Rongeur


Thin footplate / 10mm opening down

	Bite Width	Length of Shaft		
		6" / 15 cm	7" / 18 cm	8" / 21 cm
90° Down Bite 	1 mm	404-480	404-530	404-424
	2 mm	404-481	404-525	404-426
	3 mm	404-482	404-526	404-428
	4 mm	404-483	404-527	404-430
	5 mm	404-484	404-528	404-432
	6 mm	404-485	404-529	404-434



Sontec Kerrison 40° Rongeur

Thin footplate / 10 mm opening / down

	Bite Width	Length of Shaft	
		7" / 18 cm	8" / 21 cm
 40° Forward Down-Bite	1 mm	404-5301	404-553
	2 mm	404-5302	404-554
	3 mm	404-5303	404-556
	4 mm	404-5304	404-557
	5 mm	404-5305	404-558
	6 mm	404-5306	404-559



Kerrison 40° Rongeur

Black and gold ceramic coating
/ thin footplate-punch / up / 8"
(20.3 cm)

	Bite Width	Metal
		Black/Gold
40° Up-Angle 	1 mm	404-400BG
	2 mm	404-402BG
	3 mm	404-403BG
	4 mm	404-404BG
	5 mm	404-405BG



Kerrison Bayonet 40° Rongeur

Thin footplate / 10 mm opening / 8" (20 cm)
shaft / with & without ejection pin

		Bite Width	With Pin	Without Pin
40° Up-Angle		1 mm	NA	404-401
		2 mm	404-4052	404-409
		3 mm	404-405	404-4033
		4 mm	404-4054	404-407
		5 mm	404-4055	404-403
				



Kerrison Bayonet 40° Rongeur

Black ceramic coating / thin footplate / 10 mm opening / 8" (20 cm) shaft / with & without ejection pin

	Bite Width	With Pin	Without Pin
40° Up-Angle 	1 mm	NA	404-401B
	2 mm	404-4052B	404-409B
	3 mm	404-405B	404-4033B
	4 mm	404-4054B	404-407B
	5 mm	404-4055B	404-403B



Kerrison Bayonet Foraminotomy Rongeur

Curved up / 10" (25.4 cm) shaft

404-468 2 mm jaw with thin footplate
404-4683 3 mm jaw with thin footplate
404-4684 4 mm jaw with thin footplate



Kerrison Convertible 40° Rongeur


8" (20.3 cm) shaft / 10 mm opening / without
ejection pin / ejector pin for 3 mm sizes and
longer only

* Other tip sizes & lengths available upon request.

* Swing system for easy cleaning

* Sizes available 4-6 mm with pin

		Bite Width		
40° Up-Angle		1 mm	404-337C	without pin
		2 mm	404-338C	without pin
		3 mm	404-339C	with pin
		4 mm	404-340C	with pin
		5 mm	404-341C	with pin
		6 mm	404-342C	with pin



Kerrison 40° Rongeur Titanium

7" (18 cm)

2603-541 3 mm



2603-541

Sontec Kerrison 40° Rongeur

10 mm opening / up / with ejection pin

		Bite Width	Length of Shaft	
40° Up-Angle			7" / 18 cm	8" / 21 cm
		2 mm	404-540EJ	404-402EJ
		3 mm	404-541EJ	404-404EJ
		4 mm	404-542EJ	404-406EJ
		5 mm	404-543EJ	404-408EJ



Kerrison Punch (Occular-plasti) 90° Rongeur

4" (10 cm) shaft

16-2100	1.0 mm	wide up
16-2101	1.5 mm	wide up
16-2102	2.0 mm	wide up
16-2106	2.5 mm	wide up
16-2103	3.0 mm	wide up
16-2104	3.0 mm	wide down
16-2105	4.0 mm	wide down
16-2107	4.0 mm	wide up

Black ceramic coating / 4" (10 cm) shaft

16-2101B	1.5 mm	wide up
16-2102B	2 mm	wide up



16-2101



16-2103



16-2105



16-2107



16-2102



16-2104



16-2106



16-2102



16-2101B

Kerrison Punch (Ocular-plasti) 40° Rongeur

4" (10 cm) shaft

16-2096	1.0 mm wide
16-2097	1.5 mm wide
16-2098	2.0 mm wide
16-2099	3.0 mm wide


Black ceramic coating / 4" (10 cm) shaft

16-2096B	1.0 mm wide
16-2097B	1.5 mm wide
16-2098B	2.0 mm wide


	Bite Width	
	1 mm	16-2096
	1.5 mm	16-2097
	2 mm	16-2098
	3 mm	16-2099



16-2098



16-2096




16-2096B

Sontec Kerrison Rongeur

40° / easy clean / 8" (20.3 cm)

404-400E 1 mm
 404-402E 2 mm with ejector pin
 404-404E 3 mm with ejector pin
 404-406E 4 mm with ejector pin
 404-408E 5 mm with ejector pin

	Bite Width	With Pin	Without Pin
 40° Up-Angle	1 mm	NA	404-400E
	2 mm	404-402E	NA
	3 mm	404-404E	NA
	4 mm	404-406E	NA
	5 mm	404-408E	NA



404-402E

Rongeurs Tray

Anodized aluminum tray / 9" x 13" x 6"

18-3562



18-3562



Care and Maintenance for the Universal Rotating 360° Kerrison Rongeur Shafts & Handles


The Universal Rotating 360° Kerrison Rongeur shafts are comprised of four parts:

A: outside shaft with frontal jaw blade

B: Inner shaft with threading at the tip

C: Back side blade/ jaw with screw port at the end


D: spring



Slide inner shaft (B) through outer shaft (A).
Then screw jaw (C) into inner shaft (B).



Screw inner shaft (B) and jaw (C) completely together (only
after Inner shaft (B) and outer shaft (A) are together.)



To put the jaws together, fit top jaw into the grooves
of the outside shaft (A) and slide into place.


Use the grooves on the outside shaft (A) and the back
side jaw/ blade (C) to fit the two parts into position.



↑ This is the point where the jaw/ blade (C) and inner shaft (B)
are screwed together. You must make sure the two parts are
fully screwed together for the Rongeur to work properly.


INSERTING SHAFT

RELEASING SHAFT



IMPORTANT:

WHEN INSERTING NEW SHAFT INTO NEW HANDLE, YOU HAVE TO MAKE SURE THAT HANDLE IS IN **CLOSED** POSITION AND MOVING HANDLE FREE TO MOVE. OTHERWISE, THE NEW FORCE-OPEN MECHANISM WILL NOT FUNCTION!




Simultaneously with pressing the button pull the lower shaft forward to detach it from the handle.

Click Kerrison punches are suitable for all kinds of shafts: Straight, Bayonet, Curved up and down. You do not need to keep a large stock of Kerrisons, only the shafts are necessary.

This new kind of Kerrison's approach, besides significant financial savings simplifies the entire working process. Your productivity will not be affected by a very short training period. The costs of migration to this new product are very

low and the positive results are almost immediately visible.

For detailed information regarding training and available additional instruction see the contacts on the back cover of this brochure.



Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS

Sontec Instruments, Inc.
7248 South Tucson Way
Centennial, Colorado 80112 USA
Toll-free / 800.821.7496
Local / 303.790.9411
Fax / 303.792.2606

SontecInstruments.com