



 Sontec Instruments®

M.I.S. Micro Neurology

Micro Neuro 9mm Ø Large Gauge with 11.0mm Gold Knurled Grip Bayonet Instrumentation Set

Individual items sold separately

Tubular tapered shaft 2.7mm - 2.0mm Ø at tip / 360° rotatable shaft / squeeze handle / 10 3/4" (27.3 cm)

2900-0011	needle holder	1.0 mm tip	straight
2900-0031	needle holder	1.0 mm tip	curved
2900-0061	biopsy forceps	2.0 mm tip	straight
2900-0081	biopsy forceps	4.0 mm tip	straight
2900-0101	biopsy forceps	2.0 mm tip	60°
2900-0111	grasping forceps	0.5 mm tip	straight
2900-0131	forceps	0.5 mm tip	curved
2900-0161	scissors	pointed	straight
2900-0181	scissors	blunt	straight
2900-0201	scissors	pointed	curved
2900-0221	scissors	blunt	curved
2900-0241	scissors	blunt	45°



2900-0081



2900-0011 2900-0031 2900-0101 2900-0131 2900-0161 2900-0101 2900-0201 2900-0221 2900-0241

Micro Neuro Original Lightweight 5mm Ø Gauge with 9.5mm Ø Golfball Grip Bayonet Instrumentation Set

Individual items sold separately

Tubular tapered shaft 2.7mm - 2.0mm Ø at tip / 360° rotatable shaft / squeeze handle / 10" (25.4 cm)

2900-001	needle holder	1.0 mm tip	straight
2900-003	needle holder	1.0 mm tip	curved
2900-006	biopsy forceps	2.0 mm tip	straight
2900-008	biopsy forceps	4.0 mm tip	straight
2900-010	biopsy forceps	2.0 mm tip	60°
2900-011	grasping forceps	0.5 mm tip	straight
2900-013	forceps	0.5 mm tip	curved
2900-016	scissors	pointed	straight
2900-018	scissors	blunt	straight
2900-020	scissors	pointed	curved
2900-022	scissors	blunt	curved
2900-024	scissors	blunt <td>45°</td>	45°



2900-003



2900-001

2900-003

2900-010

2900-013

2900-016

2900-018

2900-020

2900-022

2900-024


Heartline Micro Neuro Bayonet Large Gauge 9mm Ø with 11mm Ø Gold Knurled Grip Needle Holder

Shaft tapered from 2.7mm Ø to 2mm at tip with 360° rotation / squeeze handle / 13" (33 cm)

2900-002	without lock
2900-004	with lock



2900-002



Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS

Sontec Instruments, Inc.
7248 South Tucson Way
Centennial, Colorado 80112 USA
Toll-free / 800.821.7496
Local / 303.790.9411
Fax / 303.792.2606

SontecInstruments.com