

SONTEC LAPAROSCOPIC BIPOLAR

FINEST QUALITY IN THE INDUSTRY / AGGRESSIVE PRICING



Claw Forceps 2695-212



Micro-Tip Curved Forceps
2695-204



Delta Grasper
2695-206



Large Grasping
Forceps 2695-202



Metzenbaum Curved
Scissors 2695-210



Maryland Curved
Dissector 2695-200



Delta Grasping Curved
Forceps 2695-208



BIPOLAR POWER GRIP & BIPOLAR POWER CUT COAGULATION INSTRUMENTATION

Universal and rotatable handpiece offers a variety
of tip options sold individually or as a set

1. Universal Handle Rotatable 360° / 2695-150
2. Maryland Curved Dissector / 2695-200
3 mm Ø / 34 cm shaft
3. Obturator Outer Shaft / 2695-110
4. 5 m Bipolar Cord Banana Plug / 2695-100

 Sontec Instruments®

A TRADITION OF FINE SURGICAL INSTRUMENTS